

User Guide

Mesh Side Rail Pads – 1275/PS & 1275/PS/EXT

SIDHIL

Introduction

Thank you for purchasing a Sidhil Mesh Side Rail Pad. The Mesh Side Rail Pad is designed to provide side rail padding whilst still allowing the occupant of the bed to be able to see through the side rail and therefore reducing the feeling of restricted space. It allows the side rails to be lowered without removing the pad.

Specifications



1275/PS		1275/PS/EXT
<ul style="list-style-type: none">Sidhil Bradshaw bed rangeSidhil Freedom II bed		<ul style="list-style-type: none">Sidhil Bradshaw bed range + platform extensionSidhil Freedom II bed + platform extension
Fabric/Mesh Colour:	Cream	
Fitting method:	Velcro	
Length:	1980mm (78")	2110mm (83")
Height:	300mm (12")	300mm (12")
Net weight per pair:	2.2kg (4.9lb)	2.3kg (5.1lb)

Fitting



- Before the pads are fitted a risk assessment must be performed on a patient by patient basis to ensure suitability (including but not limited to: risk of entrapment, asphyxiation, insecure fitting and non-compatibility).**
- Side rails pose a potential entrapment hazard, refer to side rail safety notice, SR/PSN 01/01/1, supplied with the bed.**

Ensure the Velcro fastenings are located to the outside of the bed away from the bed occupant. Ensure the Velcro fastenings are fastened along the entire length of the side rail.

The side rail pads may be removed for additional ease of raising or lowering the side rails bars if desired.

Cleaning

Wipe down with a disposable cleaning cloth, moistened with a mild detergent and diluted in hot water (60°C). Extra care should be taken around areas where excess dirt or dust may gather. Rinse down with clean water and dry off with a paper towel. Always ensure that the pad is dry before re-using.

It is not recommended to use any bleach or surface cleaner that may mark the surface.

DO ensure the side rail pads are cleaned on a regular basis

Warranty

Sidhil Limited standard terms and conditions apply to all purchases; a copy is available from the address below. Your statutory rights are not affected.

Note: Although the Mesh Side Rail Pad is designed for use with Sidhil Bradshaw & Freedom II beds, subject to assessment of suitability it may be possible to use the product with other bed and side rail combinations.



CONTACT INFORMATION

Tel: 01422 233000 Fax: 01422 233010

Email: sales@sidhil.com www.sidhil.com

Sidhil Business Park, Holmfield, Halifax, HX2 9TN

A member of the Siddall & Hilton Ltd. Group of Companies



Certificate No. FM14550

INSTRUC-1275-PS rev1 08/03/12